



Reference: PTXBC003516

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TBC1 domain family member 17		
For research use only.			

GENE DESCRIPTION	TBC1 domain family, member 17
GENE ID	79735
ORIGIN SPECIES	Human
ACCESSIONS	BC003516
VECTOR	pDONR223
NCBI SYMBOL	TBC1D17
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]