

RPS25-ribosomal protein S25 Gene

Reference: PTXBC003537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S25; 40S ribosomal protein S25		
For research use only.			

GENE DESCRIPTION	ribosomal protein S25
GENE ID	6230
ORIGIN SPECIES	Human
ACCESSIONS	BC003537
VECTOR	pDONR223
NCBI SYMBOL	RPS25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgcctaaggacgacaagaagaaggacgtggaagtcggccaagaagacaagaccagtgaacaaatccggggcgaaggccaaaagaagaagtggtccaaaggcaagtcgggacaagctcaataacttagtctgttgacaagctacctatgataaactctgta
agggaagtcccaactataaactataacccagctgtgtctctgagagactgaagattcaggctccctggccaggcgagcccttcaggagctccttagtaaggacttatcaaaactgggttcaagcacagagctcaagtaatttacaccagaataaccaagggtggagatctccagctg
ctggtgaagatgcatga