



Reference: PTXBC003552

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CNX; IP90; P90; major histocompatibility complex class I antigen-binding protein p88		
For research use only.			

GENE DESCRIPTION	calnexin
GENE ID	821
ORIGIN SPECIES	Human
ACCESSIONS	BC003552
VECTOR	pDONR223
NCBI SYMBOL	CANX
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]