

UBE2E3-ubiquitin-conjugating enzyme E2E 3 (UBC4/5 homolog, yeast) Gene

Reference: PTXBC003554

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UBCH9; UbcM2; ubiquitin-conjugating enzyme E2 E3; E2 ubiquitin-conjugating enzyme E3; ubiquitin carrier protein E3; ubiquitin conjugating enzyme E2E 3; ubiquitin-conjugating enzyme E2-23 kDa; ubiquitin-conjugating enzyme E2E 3 (UBC4/5 homolog, yeast); ubi		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2E 3 (UBC4/5 homolog, yeast)
GENE ID	10477
ORIGIN SPECIES	Human
ACCESSIONS	BC003554
VECTOR	pDONR223
NCBI SYMBOL	UBE2E3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccagtcataggcaaggtccgatgatgagagccccagcaccagcagtgccagttcagatgcggaccagcgagaccagccgtccagagcctgaagaacaagaggaaagaaaccttctgccaccagcagaagaaaaacacaaactctctagcaaaaccactgctaagtta
tccactagtctaaaagaattcagaaggagctagctgaataacccttgatcctcctctaattgcagtgctggcctaaaggagataacattatgaatggagatcaactatactggccaccgggttctgtatatgaaggtggtgtgtttttctggatatacattttcatcagattatccatttaa
gccaccaaggttactttccgaccagaatctatcactgcaacatcaacagtcaggagtcattctgttgacatccttaagacaactggagtccttgaactattcaaggttttgcgtctattgttccctttgacagactgcaaccctgcggatcctctggttgaagcatagccactc
agtatttgaccaacagagcagaacacgacagatagccagacagtgaccaagagatacgaacataa