



Reference: PTXBC003555

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOC2L nucleolar associated transcriptional repressor; NET15; NET7; PPP1R112; nucleolar complex protein 2 homolog; NOC2-like protein; novel INHAT (inhibitor of histone acetyltransferase) repressor; novel INHAT repressor; nucleolar complex associated 2 homo		

For research use only.

GENE DESCRIPTION	nucleolar complex associated 2 homolog (S. cerevisiae)
GENE ID	26155
ORIGIN SPECIES	Human
ACCESSIONS	BC003555
VECTOR	pDONR223
NCBI SYMBOL	NOC2L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]