

AP1S1-adaptor-related protein complex 1, sigma 1 subunit Gene

Reference: PTXBC003561

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AP19; CLAPS1; EKV3; MEDNIK; SIGMA1A; AP-1 complex subunit sigma-1A; HA1 19 kDa subunit; adapter-related protein complex 1 sigma-1A subunit; adaptor protein complex AP-1 subunit sigma-1A; adaptor-related protein complex 1 subunit sigma-1A; clathrin assembl		

For research use only.

GENE DESCRIPTION	adaptor-related protein complex 1, sigma 1 subunit
GENE ID	1174
ORIGIN SPECIES	Human
ACCESSIONS	BC003561
VECTOR	pDONR223
NCBI SYMBOL	AP1S1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcggttcatgctattatcagccggcagggaactgcggctgcaaaaatggtaacctggccactcggacaagggaagagatggcgagctcatgcaggttgctcggctcgaagcccaagatgtgcagcttctggagtgaggacacaaagttgtctataag
agatatgccagcctctacttctgctgcgcacatcaggcccaagacaatgagctcatcacactggagctgatccaccgatacgtggagctcttagacaatacttggcagtggtgagctggacatcatctcaactttgagaaggcctactcctgagatgtttgatggggggga
tgtccaggacacctccactttcccttttccactga