



Reference: PTXBC003575

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GLCNAC-TI; GLCT1; GLYT1; GNT-1; GNT-I; GnTI; MGAT; alpha-1,3-mannosyl-glycoprotein 2-beta-N-acetylglucosaminyltransferase; N-glycosyl-oligosaccharide-glycoprotein N-acetylglucosaminyltransferase I; glcNAc-T I; mannosyl (alpha-1,3)-glycoprotein beta-1,2-N		

For research use only.

GENE DESCRIPTION	mannosyl (alpha-1,3-)-glycoprotein beta-1,2-N-acetylglucosaminyltransferase
GENE ID	4245
ORIGIN SPECIES	Human
ACCESSIONS	BC003575
VECTOR	pDONR223
NCBI SYMBOL	MGAT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]