



Reference: PTXBC003578

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FIGLER1; SALM4; leucine-rich repeat and fibronectin type-III domain-containing protein 3; fibronectin type III, immunoglobulin and leucine rich repeat domains 1; synaptic adhesion-like molecule 4; leucine rich repeat and fibronectin type III domain contain		

For research use only.

GENE DESCRIPTION	leucine rich repeat and fibronectin type III domain containing 3
GENE ID	79414
ORIGIN SPECIES	Human
ACCESSIONS	BC003578
VECTOR	pDONR223
NCBI SYMBOL	LRFN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]