

C17orf62-chromosome 17 open reading frame 62 Gene

Reference: PTXBC003595

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C17orf62; chromosome 17 open reading frame 62		
For research use only.			

GENE DESCRIPTION	chromosome 17 open reading frame 62
GENE ID	79415
ORIGIN SPECIES	Human
ACCESSIONS	BC003595
VECTOR	pDONR223
NCBI SYMBOL	C17orf62
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtacctcagggtggagacccgcaccagctcccgctccatctgaagaggggtccaggcatccggctcctgtggttggaattgtgcattggcctggctgctgcctactacagcggagatagcctgggctggaagctcttctacgtcacaggctgcctgtttgtggctgtgcag
aacttggaggactgggaggaagccatcttcacaagagcacagggaaggtgttttgaagacgttcagcctctacaagaagctgctgactctttcagagctggccacgaccaggtggtgtcctgctccatgatgtccgtgatgtgagcgtggaggagagaaggtccggtacttcggga
aaggctacatggtgtgtcctcggttcgcacgggtctccacccctcagcagagtgcatgcatgggcccaccgcagtgatgtggaagccatcgaagctcaccagcttcttgagctgactgccttgagagcccccacagagctgtctcagagcagcagtgaggccggt
gacctgcaagccagagctga