

LMO4-LIM domain only 4 Gene

Reference: PTXBC003600

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LIM domain transcription factor LMO4; LIM domain only protein 4; LIM-only 4 protein; LMO-4; breast tumor autoantigen; LIM domain only 4		

For research use only.

GENE DESCRIPTION	LIM domain only 4
GENE ID	8543
ORIGIN SPECIES	Human
ACCESSIONS	BC003600
VECTOR	pDONR223
NCBI SYMBOL	LMO4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaatccggcagcagctcgcagccgccccggtagcgccggctcctctcctggaagcgggtgcgaggtcgggggcaagattcggaccgctttctctatgccatggacagctattggcacagccgggtgcctcaagtgtcctctgctccaggcgcagctggggcagac
atcggcagctcctgttacacaaaagtgcatgatccttgcagaaatgactacattaggttatttggaaatagcgggtccttgcagcgttgcggacagtcgattcctgcgagtgaaactgtcatgagggcgcaaggcaatgtgtatcatcttaagtgtttacatgctctacctgccggaatcgc
ctggtcccgaggagatcggtttcactacatcaatggcagttattttgtgaacatgatagacctacagctctcatcaatggccatttgaattcacttcagagcaatccactactgccagaccagaaggtctgctaa