

GEMIN8-gem (nuclear organelle) associated protein 8 Gene

Reference: PTXBC003607

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FAM51A1; gem-associated protein 8; family with sequence similarity 51, member A1; gemin-8; gem nuclear organelle associated protein 8		

For research use only.

GENE DESCRIPTION	gem (nuclear organelle) associated protein 8
GENE ID	54960
ORIGIN SPECIES	Human
ACCESSIONS	BC003607
VECTOR	pDONR223
NCBI SYMBOL	GEMIN8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcggtaaagcatcaacatgaaagctaccagccttggtattctatccggtatatgcaagatactggcaacattatcatcaagcaatggcttgatgcaagccatcacaatgcctacaggagggcgtggaatcctgtttcaattctccatgggtacttactctctgcttctcc
caaagctcttacgataatgaggtgctgctcagtccttctatgaccatcatgtggcctggcaggactacccctgcagttcttcacatttcagaagatctgggcagcatccagcttacagcagtaggagccagccatccacaaaagaagaccaagctttgtccaaagaggagagatggag
actgagtcagatgcagaggtagaatgtgacctgagcaaatatggaatcactgaagagctccgccagcttctgagagaccgagagcatagagaagaacgacggcgcagcagcagctggatgcagagcgctggagagctatgtgaacgctgaccacgacctgtactgcaacacc
cgccggtcggtagaagcccaactgagagcctggtgagcggcgccaggccgagatgaagcgtttgtacggggacagtgctgccaagatccaagccatggaggccggtgagctgagctttgacaagcactgtgaccgaaagcagcccaagtactggccggtcatccctgaa
gttctga