

FAM103A1-family with sequence similarity 103, member A1 Gene

Reference: PTXBC003627

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM103A1; C15orf18; HsT19360; RAM; RNMT-activating mini protein; family with sequence similarity 103 member A1		

For research use only.

GENE DESCRIPTION	family with sequence similarity 103, member A1
GENE ID	83640
ORIGIN SPECIES	Human
ACCESSIONS	BC003627
VECTOR	pDONR223
NCBI SYMBOL	FAM103A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactgacactgccgaagctgttccaagtttgaagagatgttctctagtagattcacagaaatgacaaggagatcaggaatacctgaacgccctcctgagtcctccaattgttgaggatggaatagcagagctggtgggaaccaagaacagaggcaatcggttgcagacaa
cagacagttcagaggcagggaacaacagatgggggtggccaagtgaacatcgatccaatcagtgccatggacatcctgggtaacaactaccgcaacacagacaagaaccttactatccccagcaatattgacattatggttacaaccagcggcctccttacggttactactga