

C16orf53-chromosome 16 open reading frame 53 Gene

Reference: PTXBC003640

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C16orf53; GAS; PA1; PAXIP1-associated glutamate-rich protein 1; PAXIP1-associated protein 1; PTIP-associated 1 protein; PTIP-associated protein 1; glutamate-rich coactivator associated with SRC1; glutamate-rich coactivator interacting with SRC1; glutamate		

For research use only.

GENE DESCRIPTION	chromosome 16 open reading frame 53
GENE ID	79447
ORIGIN SPECIES	Human
ACCESSIONS	BC003640
VECTOR	pDONR223
NCBI SYMBOL	C16orf53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcccttctcggggccatggagacactgcggccagtcggggcgccctctgtctgaagaaggggaagtgcctccggcctccagctctggcgtggagataccggagccctctgcctcggccggttaaggccgagggagggggaaggaggccgagaggagaccga
gcctgaggggtccggggcgaggagggcgaggagagaagtcgccagcgtgggggagaagagcctgccgaggaggactccgaggactggctgcctgcctgcagcgacgagggagtgagctgcctgcctggatggcgagccctggatgccccgcctccgaaatccagcggct
ctatgaactgctggctgcccacgggtactctggagctgcaagccgagatcctgccccgccgcccctccacgcccggaggcccgagcgaagaggagagatccgatgaggagccggagggccaaagaggaagaggaagaaaaccacacatgcccacggaatttgatttgatgatgag
ccagtgcaccaaaggactccctgattgaccggagacgcacccaggaagctcagccggagccagaacgggaggcccgctggacaaggctgctgtcgacatgaagagacacaagagctggaggagcagatccttcgtaccggaggagacctctcagcctgactcggag
gacccagcccccgagccccccacfcgatcctccggagtagtctcttccctcggcagcggaataactga