

WDR25-WD repeat domain 25 Gene

Reference: PTXBC003641

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf67; WD repeat-containing protein 25; WD repeat domain 25		
For research use only.			

GENE DESCRIPTION	WD repeat domain 25
GENE ID	79446
ORIGIN SPECIES	Human
ACCESSIONS	BC003641
VECTOR	pDONR223
NCBI SYMBOL	WDR25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagctctgtgctggtgctgtgctgttcaagggtgttggagaggtatggaacgccgtggactccgggcaactgcagacctactccctgcacacagaggcagtcggggccggcggtgggtccctgtggccggcgcatcctcagtggtgctttgacttcgcgtgcac
ctaacagaccttgaacaggaaccagctatttagtgcgaagtgaacttgaatcactaccttgaattccatccaaaagaccacaacatcttttatgtggaggttcagctctgaaatgaaagcttggatataaggaactggcaaggtgatgagaagctacaaggcgacctccagcaga
ccttggacatcctgttctccgggaaggctccgagttcctgagcagcacagacgttccacccgggactcagctgaccgcaccattattgcctgggatttccggacctctgccaaatctccaaccagattttccacgagaggttcacctgccccagcctgccttgcacccgagagagccc
gtgttctgtgcacagaccaatggcaactacgtggcccttttctcactgtgtggccctaccggatgagcagacggcgcgctatgaagggcacaaggtggagggctactcagtggtgctgcgagtgctccccagcgcggtgacttggtgacgggcagcgccgatggccggtcctga
tgtacagcttccgcacagccagccgagcatgcacactgcaggggcacacacagggcctgtgtggcaccacctaccaccccgctgtgccctccgtctcgcacctgctcctggggaggggacatgaagatctggcactga