

FABP4-fatty acid binding protein 4, adipocyte Gene

Reference: PTXBC003672

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A-FABP; AFABP; ALBP; HEL-S-104; aP2; fatty acid-binding protein, adipocyte; adipocyte fatty acid binding protein; adipocyte lipid-binding protein; adipocyte-type fatty acid-binding protein; epididymis secretory protein Li 104; fatty acid binding protein 4		

For research use only.

GENE DESCRIPTION	fatty acid binding protein 4, adipocyte
GENE ID	2167
ORIGIN SPECIES	Human
ACCESSIONS	BC003672
VECTOR	pDONR223
NCBI SYMBOL	FABP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgatgcttttaggtacctggaaactgtctccagtgaactttgatgattatgaagaagtaggagtggtgcttgcaccaggaaagtggctggcatggccaaacctaactgatcatcagtgatggtgggatgtgatcaccattaaatctgaaagtacctttaaatactgaga
tttcctcatactgggcccaggaaatttgacgaagtcactgcagatgacaggaaagtcaagagcaccataacttagatgggggtgtctgtacatgtgcagaaatgggatggaaatcaaccaccataaagagaaaacgagagatgataaactggtgtggaatgcgtcatgaaagcgt
cactccagagagtttatgagagacataa