

PRM1-protamine 1 Gene

Reference: PTXBC003673

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CT94.1; sperm protamine P1; cancer/testis antigen family 94, member 1; cysteine-rich protamine; testicular tissue protein Li 91; testis specific protamine 1; protamine 1		

For research use only.

GENE DESCRIPTION	protamine 1
GENE ID	5619
ORIGIN SPECIES	Human
ACCESSIONS	BC003673
VECTOR	pDONR223
NCBI SYMBOL	PRM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaggtacagatgctgtcagccagagccggagcagatattaccgacagacaagaagtcgacagcaaggaggcgagctgccagacacggaggagagccatgaggtgctgccgccccaggtacagaccgagatgtagaagacactaa