

NDUFA2-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 2, 8kDa Gene

Reference: PTXBC003674

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD14; CIB8; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 2; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 2, 8kDa; NADH-ubiquinone oxidoreductase subunit CI-B8; complex I B8 subunit; NADH:ubiquinone oxidoreductase subunit A2		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 2, 8kDa
GENE ID	4695
ORIGIN SPECIES	Human
ACCESSIONS	BC003674
VECTOR	pDONR223
NCBI SYMBOL	NDUFA2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgccgcagcaagtcgaggagtcggggcaagctgggcctgctgagattcgatccactatgtcagcgctgccccgcagccaggcgctcagggacttcattgagaaacgctacgtggagctgaagaaggcggaatcccgaacctaccatcctaaccggaatgctcc
gatgtgcagcccaagctctggcccgctacgcatttggccaagagacgaatgccccttgaacaacttcagtgtgatcaggtaccagagccctggagaacgttctaagtggtaaagcctga