

ATP5I-ATP synthase, H⁺ transporting, mitochondrial F0 complex, subunit E Gene

Reference: PTXBC003679

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP5K; ATP synthase subunit e, mitochondrial; ATP synthase e chain, mitochondrial; ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit E; ATPase subunit e; F1F0-ATP synthase, murine e subunit; ATP synthase, H ⁺ transporting, mitochondrial Fo c		

For research use only.

GENE DESCRIPTION	ATP synthase, H ⁺ transporting, mitochondrial F0 complex, subunit E
GENE ID	521
ORIGIN SPECIES	Human
ACCESSIONS	BC003679
VECTOR	pDONR223
NCBI SYMBOL	ATP5I
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgccaccgggtgcaggtctctccgtcatcaagctcggccgtactccgctgttctcgggtgtggcctacggagccacgcgtacaattacctaaaaacctcgggcagaagaggagaggagtagcagcagaagagaagaagacaggatgaactgaacggattgccaga
gaattggcagaagatgacagcatattaaagtga