



Reference: PTXBC004141

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGAT; CCDS3; glycine amidinotransferase, mitochondrial; glycine amidinotransferase (L-arginine:glycine amidinotransferase); testicular secretory protein Li 19; transaminidase; glycine amidinotransferase		

For research use only.

GENE DESCRIPTION	glycine amidinotransferase (L-arginine:glycine amidinotransferase)
GENE ID	2628
ORIGIN SPECIES	Human
ACCESSIONS	BC004141
VECTOR	pDONR223
NCBI SYMBOL	GATM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]