



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KRTAP4-12-keratin associated protein 4-12 Gene

Reference: PTXBC004180

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KAP4.12; KRTAP4.12; keratin-associated protein 4-12; keratin-associated protein 4.12; ultrahigh sulfur keratin-associated protein 4.12; keratin associated protein 4-12		

For research use only.

GENE DESCRIPTION	keratin associated protein 4-12
GENE ID	83755
ORIGIN SPECIES	Human
ACCESSIONS	BC004180
VECTOR	pDONR223
NCBI SYMBOL	KRTAP4-12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]