



Reference: PTXBC004190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CBARA1; CALC; EFHA3; MPXPS; calcium uptake protein 1, mitochondrial; ara CALC; atopy-related autoantigen CALC; calcium-binding atopy-related autoantigen 1; mitochondrial calcium uptake 1		

For research use only.

GENE DESCRIPTION	calcium binding atopy-related autoantigen 1
GENE ID	10367
ORIGIN SPECIES	Human
ACCESSIONS	BC004190
VECTOR	pDONR223
NCBI SYMBOL	CBARA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr