

DCTN5-dynactin 5 (p25) Gene

Reference: PTXBC004191

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	dynactin subunit 5; dynactin 4; dynactin 5 (p25); dynactin subunit p25		
For research use only.			

GENE DESCRIPTION	dynactin 5 (p25)
GENE ID	84516
ORIGIN SPECIES	Human
ACCESSIONS	BC004191
VECTOR	pDONR223
NCBI SYMBOL	DCTN5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagttggcgagctgctctacaacaagtctgagtagcatcgagagggcatctgggaacaagtcagtcgccagtcagtggtgtgtggaagccagaacatggtctcaatggcaagaccattgtgatgaatgactgtattatccgaggggatctggcaaatgtaagagttggacgtcattgtgtgtgaaaagtcgtagtgataaggccaccattcaagaagttcagcaaaaggtgtgcatcttctttacataattggagaccatgctttatggagagattgtgtgtgcaacgcagcacagattggttcctatgttcatgttgggaagaactgtgtgattgggcgccgatgtgtgtgaaagacgtctgcaaaattcttgacaacacagtattacctccagaactgtgtgtccaccattcactgtcttctcaggctgccaggactcttctcagggagctcccggagtgacactcaggagctgatgattgacgtcaccaagagctactaccagaagttttgccctgacgcaagtcctag