

HNRNPUL1-heterogeneous nuclear ribonucleoprotein U-like 1 Gene

Reference: PTXBC004242

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E1B-AP5; E1BAP5; HNRNPUL1; heterogeneous nuclear ribonucleoprotein U-like protein 1; E1B 55kDa associated protein 5; E1B-55 kDa-associated protein 5; adenovirus early region 1B-associated protein 5; heterogeneous nuclear ribonucleoprotein U like 1		

For research use only.

GENE DESCRIPTION	heterogeneous nuclear ribonucleoprotein U-like 1
GENE ID	11100
ORIGIN SPECIES	Human
ACCESSIONS	BC004242
VECTOR	pDONR223
NCBI SYMBOL	HNRNPUL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgtgcgcgtctgaaggtgaacgaacttcgcgaggagctgcagcgcgcggcctggacactcaggcctcaaggccgagcttgctgagcggctgcaggcggcgttgaggccgaggagcctgacgacgagcgggagctcagcgcgacgacgaacggggcgacccg
ggcacatcaacgaggagccgaggttacagcagccacccctacaaccaggaggttacagccagggctacacagcccaccgctccacctccaccaccactgcctacaactatgggagctacggcggttacacccggcccccctatacccccaccgccccccaccgcacaga
cctaccctcagccagctataa