



Reference: PTXBC004249

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	OIP-1; OIP1; TRIP-6; TRIP6i2; ZRP-1; thyroid receptor-interacting protein 6; OPA-interacting protein 1; TR-interacting protein 6; thyroid hormone receptor interacting protein 6; zyxin related protein 1; thyroid hormone receptor interactor 6		

For research use only.

GENE DESCRIPTION	thyroid hormone receptor interactor 6
GENE ID	7205
ORIGIN SPECIES	Human
ACCESSIONS	BC004249
VECTOR	pDONR223
NCBI SYMBOL	TRIP6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]