

CBX2-chromobox homolog 2 (Pc class homolog, Drosophila) Gene

Reference: PTXBC004252

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDCA6; M33; SRXY5; chromobox protein homolog 2; Pc class homolog; cell division cycle associated 6; chromobox homolog 2 (Pc class homolog, Drosophila); modifier 3; chromobox 2		

For research use only.

GENE DESCRIPTION	chromobox homolog 2 (Pc class homolog, Drosophila)
GENE ID	84733
ORIGIN SPECIES	Human
ACCESSIONS	BC004252
VECTOR	pDONR223
NCBI SYMBOL	CBX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagctgagcagctggcgagcaggtcttcgccgagtgatcctgagcaagcggtccgcaaggcgaagctggagtacgtggtcaagtggcgcggtctggtctccaataacagctgggagccggaggagaacatcctggaccgaggtgctcctggtccca
gaagaaggaaatgagaaggaggtgcagaaccggaagagaggcaagagccgagagccgccaaggaaagctcactgcatgtctctgcaagccggcgctccaagctcaaggtgggtggctgctggtatgtagacccacacccagcacccttggcgtaggggca
ggcagagggaggtttggggccctcaggaagggtgagcacttctgccaacagctgtctcctctactcggaaaacagagagcccttctctctctcagcttctgctgccagggggcccgccgctgagagttctcccccctccaggggtctctcagcctgtctgc
acgcctctctgctgggtggcagggtcaactgctgcagacaagcactcttccctccaggggtccttggggagcggaaaggaacaggaagcatgctacagtag