

FAM164C-family with sequence similarity 164, member C Gene

Reference: PTXBC004259

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM164C; FAM164C; C14orf140; zinc finger C2HC domain-containing protein 1C; family with sequence similarity 164, member C; zinc finger C2HC-type containing 1C		

For research use only.

GENE DESCRIPTION	family with sequence similarity 164, member C
GENE ID	79696
ORIGIN SPECIES	Human
ACCESSIONS	BC004259
VECTOR	pDONR223
NCBI SYMBOL	FAM164C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtctccagcggttgccgtcacatctgctgtggcggttatgtcccacataatacaacggaagctccaggccccactcagccaagcaagactttacgaacaaggtgactcttccagcagtcctgaaggggcacctgaggacaattccagaagcagcttttgagcaaca
aagagttgattctggataaagtctatactcaccacaaatgggaacccaacaaaagccggagctactctatccccactgtactggaaatcagccagcaagatccagaagtgatgccaggccaaggaatggtttgttttactgtcaggccctcaatcctggtatccaaagccaata
accaggactttatcccccttacaagaacagagttggagtgaccggcggttccatgaaacccatggtccacaggaagctgtgagtagcaggtgaggtggcactgatggggaccataatgtctaccagggccctgagccgagagagttttcatcaggaaactttggtgtgaggaa
ccagggcaacttttctgtgtgtgtactgtctgtgtgcccacgcagggcgagaaggccgtggcacaactttgacaggacggagtggtgtgcagatccgaagactagaagctgcagggggagagcttagaggaggaaatccgaagaaagcagattctctgaggggaaagctgaaagagac
agaggaggaaactcagaaggatccagacgcaaaaggacagggccaaggaaaatgaaaacggagagctacagaaaattatctcccaggagcagagtttaaaggctga