

PRR7-proline rich 7 (synaptic) Gene

Reference: PTXBC004261

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	proline-rich protein 7; synaptic proline-rich membrane protein; proline rich 7, synaptic		
For research use only.			

GENE DESCRIPTION	proline rich 7 (synaptic)
GENE ID	80758
ORIGIN SPECIES	Human
ACCESSIONS	BC004261
VECTOR	pDONR223
NCBI SYMBOL	PRR7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgatgtccaggacacgttctcagctgttctgcccgttctggtctatctggtctctgctgttctgagcttctgcccgcgcctcaaacggcgccaggaggagcgtgcccgcgcagaaacctgcgcgccctagagctggagccctcgaac
tcgaggcgagctgtgcccggagccccccggcctggcgccgcgcagccaccacacacgtagccgctggaggcgccggctcagcgccgcgcacccacacccgcaccacacgcgtcccgacccgcgcctacgcacctgtcggfcccgcacggccctggagta
cccgcgcgaagcgggaatcggacatgtcacaaccacgtgttacgaaggcggtgctgatggcagagccgcgcgccctatagcagggtgctcagcgacacgcgcgcctctaccgcaagatcgtacgcccctctgagtcgccgcgacagcgcggagaagcaggagcagcc
gcctcccagctacaagccgcttctcgtggaccggggtacacctggcgctgcacctgcccagcgcgccctggcccgcgcgcctgcccagccctgctgctgagggcgacccgtggccgggtcttccccagctggaccgactcagagctcagcagccgcgagccctggagc
acggagcttggcgtctgcgggtctccatccctgttggaggagactacagccgtatag