

CERK-ceramide kinase Gene

Reference: PTXBC004278

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LK4; dA59H18.2; dA59H18.3; hCERK; ceramide kinase; acylsphingosine kinase; lipid kinase 4; lipid kinase LK4		

For research use only.

GENE DESCRIPTION	ceramide kinase
GENE ID	64781
ORIGIN SPECIES	Human
ACCESSIONS	BC004278
VECTOR	pDONR223
NCBI SYMBOL	CERK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaactggccctgctcaggtgatccggggaggggctccacgcagccgggaaggctgctggggcctccacacgtgctcatcacggcgaggctacgacaatccggctgggagcatgacctggcgctgtcttgggagcacggatgataagctctggaagctggca
gtgtgtaaagcactggcgaagtgtgtactgttaaatgtcaaataccaatgctttatcgacgcgaagtgtctaacacagccggcttgggggagtcaggagggaagctggccatccgtggaggagggcggtcctggactcccgaggactcctctgaggcaggccctgaagtctgta
cacgtgtccagattgtcctgtcttcttcacactgagttctctatattatgaacatctgtctctttaagccagagtagtgtaaactgcgtctcggatgtctgtctttgctcgaagccacgatggatcgctgggttctctgcagcgcgagggtcggcgaccagaggattctccggaag
gcattcctgccgcgtcccggggacccctcaattgtgtactacgtcctgtttag