



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SYT11-synaptotagmin XI Gene

Reference: PTXBC004291

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

SYT12; sytXI; synaptotagmin-11; synaptotagmin 12; synaptotagmin XI; synaptotagmin 11

For research use only.

GENE DESCRIPTION	synaptotagmin XI
GENE ID	23208
ORIGIN SPECIES	Human
ACCESSIONS	BC004291
VECTOR	pDONR223
NCBI SYMBOL	SYT11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr