

PSMG3-proteasome (prosome, macropain) assembly chaperone 3 Gene

Reference: PTXBC004308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf48; PAC3; proteasome assembly chaperone 3; PAC-3; hPAC3; proteasome (prosome, macropain) assembly chaperone 3; proteasome assembling chaperone 3		

For research use only.

GENE DESCRIPTION	proteasome (prosome, macropain) assembly chaperone 3
GENE ID	84262
ORIGIN SPECIES	Human
ACCESSIONS	BC004308
VECTOR	pDONR223
NCBI SYMBOL	PSMG3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagacacgccgttggtgatacgaagcagaagacggaggtgtgtgcgggggtcccccacccaggtgtgtgtacggccttcagcagtcacatcctgtgtgtgtgacccagtttgggaagatgggcaccctgtctccctggagcccagcagcgtggccagtgcagtcagcaagc
ctgtgtcaccacaaaagtcctctgtgggcaggatgagcctctcatccatgtctttgcaaaagacctgtgagcgtttgtgtctcaagaagctggaaacagagcagtcctctcgcctgtggcgtgaaggacaaaagcatggagggtgaaggcgtgaggaggtgatccgggtgtgcc
aggtgtgtgtga