

TMBIM4-transmembrane BAX inhibitor motif containing 4 Gene

Reference: PTXBC004401

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-119; GAAP; LFG4; S1R; ZPRO; protein lifeguard 4; Golgi anti-apoptotic protein; Z-protein; transmembrane BAX inhibitor motif-containing protein 4; transmembrane BAX inhibitor motif containing 4		

For research use only.

GENE DESCRIPTION	transmembrane BAX inhibitor motif containing 4
GENE ID	51643
ORIGIN SPECIES	Human
ACCESSIONS	BC004401
VECTOR	pDONR223
NCBI SYMBOL	TMBIM4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgaccccgaccccggtaccctcgctcctcgatcgaggacgacttcaactatggcagcagcgtggcctccgccaccgtgcacatccgaatgggtctcttaactacagtgaactcaacagtttttatacttga