

TMEM177-transmembrane protein 177 Gene

Reference: PTXBC004404

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 177		
For research use only.			

GENE DESCRIPTION	transmembrane protein 177
GENE ID	80775
ORIGIN SPECIES	Human
ACCESSIONS	BC004404
VECTOR	pDONR223
NCBI SYMBOL	TMEM177
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagggtccctgtggcggaccgcagcattgtgcagagacacaggacggcctctgtgggttcctgtgcaggcctgtttggagttccagtcgtaccacctcttccggatcccggtggtccaatggctctaccagtactggcctcaggccagccagctccctccacagct
gcagagcctctccaagggtgtctacaggacataggtgtcttcaggccattgtctacaagcccttaccaccttccagcctgtgagtgaggcttcccaagactccctgtgggctgtgtgggcatccctgccagttcttgggagacctagtatcaactaaccatccctgtg
gtcatacatgggcatacagtgactggcggagccagcaggcggcggctgagagcttccctgacctgtgcccgtgaagccagaagttgccttggccagggaagtggtgtacctggaaagcagttaccactgcccgtgcacgcccctgtgccccagcttgcctggcagggacctgg
gcactgggctgtgggtccaagtacacctggggtccatgcaggcccatgaatttacgggtgccttcagcttgggtggcagcagtgggcaggcttgtgtggcctacgcttctccaggattctctactcatgcccgtggagtcctgtgctggaccggcgcacggcctccctctctgcagccta
tgccctgtgggtggagtggttctatgagaagcttctgtcgggcaacctggccctgcgcagttcttgggcaagagggggagaagctgtatataccacagcgggaacatgtccccagacactgttccgaatcaaacatttaccctacaccaccgccgggactctgtgtgcagatgtgg
agggggatgctcaatccgggcccgtcctga