



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

TOLLIP-toll interacting protein Gene

Reference: PTXBC004420

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-1RAcPIP; toll-interacting protein; adapter protein; toll interacting protein		
For research use only.			

GENE DESCRIPTION	toll interacting protein
GENE ID	54472
ORIGIN SPECIES	Human
ACCESSIONS	BC004420
VECTOR	pDONR223
NCBI SYMBOL	TOLLIP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgcaccaccgctcagcactcagcgcggcggtgtatctcgttgtagctccgcaggactctctccgtacacgccacacagcagcagcggcaggctccagctggacgccagcgcgccagcagctgcagfacaggagcgacagtggccaccgtggcccgactgaacatcacg
gtgtgtacaggccaagtgtggccaagaattacggcatgaccgcgatggaccctactctccgactgcgcctgggctctacgcggtgtacgagacgccacggcacacaatggcgccaagaatccccctgtgaaataggctacactctcagcgggtccccccaggcgtgactctttctatctga
gatcttcgatgagagagactcttcctatggacgaccgcatctgctggaccacacatccatccagagctccctgaggcaggcgcaaggtggagacaagaatgtgtacacgcttgacggcgaggcgaggcgaggacgacaagaaggcgcgatgataacactctctatctctacgcgctgtccacgc
tgccatggatgatgccaccccgccggctgtgatgcccaagctgtacacagcggcgtgtggctatgtgcccaacaggatggcctgtctgtctgacccccgcgatgtgccccgtgccccgtccccccgcgcgcgtgaacgccagcccgctgtgagcgagcggctgaagc
catccagacagctcttcccccaatctgaccagagcgtgatccctctctctctgaaagcccaagcgggaacaagatgcgcgccataacctctctctcagatggggagcccgccatg