



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CNPY3-canopy 3 homolog (zebrafish) Gene

Reference: PTXBC004423

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAG4A; ERDA5; PRAT4A; TNRC5; protein canopy homolog 3; CAG repeat containing; CTG repeat protein 4a; expanded repeat-domain protein CAG/CTG 5; protein associated with TLR4; protein associated with Toll-like receptor 4A; trinucleotide repeat-containing gen		

For research use only.

GENE DESCRIPTION	canopy 3 homolog (zebrafish)
GENE ID	10695
ORIGIN SPECIES	Human
ACCESSIONS	BC004423
VECTOR	pDONR223
NCBI SYMBOL	CNPY3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr