

NINJ1-ninjurin 1 Gene

Reference: PTXBC004440

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NIN1; NINJURIN; ninjurin-1; nerve injury-induced protein-1; ninjurin 1		

For research use only.

GENE DESCRIPTION	ninjurin 1
GENE ID	4814
ORIGIN SPECIES	Human
ACCESSIONS	BC004440
VECTOR	pDONR223
NCBI SYMBOL	NINJ1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactgggaaccgaggagtagcagctcaacggcgccctgctcgggcacacccggctcccggaacgctcgcggcccgctggggctggaggcacgggccatcaacgtgaaccattacgccagcaagaagagcgcagccgagagcatgctggacatcgcgctgctgat
ggccaacgcgtccagctgaaggccgctggaacagggccccagcttcgcttctatgtgcccctgggtggtcctcatctccatctccctgtgctgcagatcggcggtgggtgctgctcatctcttgaagtacgaccttaacaacccggccaagcacgccaagctggacttctcaa
caacctggccacggcgctggtttcatcatcgtggtagtagcaacatctcatcacggccttcgggtccagaagccctgatggacatggcaccacagcagtag