

DHX37-DEAH (Asp-Glu-Ala-His) box polypeptide 37 Gene

Reference: PTXBC004463

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DDX37; DEAD/DEAH box helicase DDX37; DEAD/H (Asp-Glu-Ala-Asp/His) box polypeptide 37; DEAH (Asp-Glu-Ala-His) box polypeptide 37; DEAH box protein 37; DEAH-box helicase 37		

For research use only.

GENE DESCRIPTION	DEAH (Asp-Glu-Ala-His) box polypeptide 37
GENE ID	57647
ORIGIN SPECIES	Human
ACCESSIONS	BC004463
VECTOR	pDONR223
NCBI SYMBOL	DHX37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcttcgtggcctgctctcagggaagtgggtcaagccctgggaacctcatccatgagagctcgatcccgatgaagggtgctgccgccgtgccatctggccccgggggtga