

C10orf116-chromosome 10 open reading frame 116 Gene

Reference: PTXBC004471

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf116; AFRO; APM2; apM-2; adipogenesis regulatory factor; adipogenesis factor rich in obesity; adipose most abundant gene transcript 2 protein; adipose specific 2; adipose-specific protein 2		

For research use only.

GENE DESCRIPTION	chromosome 10 open reading frame 116
GENE ID	10974
ORIGIN SPECIES	Human
ACCESSIONS	BC004471
VECTOR	pDONR223
NCBI SYMBOL	C10orf116
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaagcaagggcttcaggacctgaagcaacaggtggaggggacccaggaagccgtgtcagcggccggagcggcagctcagcaagtggtggaccaggccacagaggcgggcagaaagccatggaccagctggccaagaccaccaggaaccatcgacaagact
gctaaccaggccctctgacaccttctctggatcgggaaaaattcggcctcctgaaatga