

CHCHD5-coiled-coil-helix-coiled-coil-helix domain containing 5 Gene

Reference: PTXBC004498

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C2orf9; MIC14; coiled-coil-helix-coiled-coil-helix domain-containing protein 5; mitochondrial intermembrane space cysteine motif protein of 14 kDa homolog; coiled-coil-helix-coiled-coil-helix domain containing 5		

For research use only.

GENE DESCRIPTION	coiled-coil-helix-coiled-coil-helix domain containing 5
GENE ID	84269
ORIGIN SPECIES	Human
ACCESSIONS	BC004498
VECTOR	pDONR223
NCBI SYMBOL	CHCHD5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcggccctagaggtcaccgctcgtactgtggccgggagctggagcagatggccagtggtggcggccaagccggaatcctggcagcgggactgtcactacctaagatgagcattgccagtgacatcctcccaccaatcatccgcagatccgccagcctgtgtca
gcctttgaggcctcgaggagtgcttcgacagaacgagcagcgtgtggcgaactgtgcagagcatatgcgcgttctcagtgccgtgagcaggtgcagccgccagctcacctgcaactgtggagcagccacttctcctctga