



Reference: PTXBC004529

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1-AGPAT2; BSCL; BSCL1; LPAAB; LPAAT-beta; 1-acyl-sn-glycerol-3-phosphate acyltransferase beta; 1-AGP acyltransferase 2; 1-AGPAT 2; 1-acylglycerol-3-phosphate O-acyltransferase 2 (lysophosphatidic acid acyltransferase, beta); lysophosphatidic acid acyltran		

For research use only.

GENE DESCRIPTION	1-acylglycerol-3-phosphate O-acyltransferase 2 (lysophosphatidic acid acyltransferase, beta)
GENE ID	10555
ORIGIN SPECIES	Human
ACCESSIONS	BC004529
VECTOR	pDONR223
NCBI SYMBOL	AGPAT2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]