

ATP6V0C-ATPase, H⁺ transporting, lysosomal 16kDa, V0 subunit c Gene

Reference: PTXBC004537

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6C; ATP6L; ATPL; VATL; VPPC; Vma3; V-type proton ATPase 16 kDa proteolipid subunit; ATPase, H ⁺ transporting, lysosomal 16kDa, V0 subunit c; H(+)-transporting two-sector ATPase, 16 kDa subunit; V-ATPase 16 kDa proteolipid subunit; vacuolar ATP synthase		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal 16kDa, V0 subunit c
GENE ID	527
ORIGIN SPECIES	Human
ACCESSIONS	BC004537
VECTOR	pDONR223
NCBI SYMBOL	ATP6V0C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccgagtcgaagagcgccccgagtagcttctgttttgcgcgtcatggcgccctcgccgccatgtcttcagcgccctggcgctgctatggcacagccaagagcggtaccggcattgcggccatgtctgtcatgcggccggagcagatcatgaagtcacatcccatgtgtca
tggctggcatcatcgccatctacggcctgggtgggtggcagtcctcatcgccaactccctgaatgacgacatcagcctctacaagagcttctccagctggcgccggcctgagcgtgggcctgagcgccctggcagccggctttgccatcgccatcgtggggacgctggcgtgcggg
caccgccagcagccccgactattcgtggcatgatcctgattctcatcttcgcgaggtgctcgcctctacgggtctcatcgtcgcctcatcctctccacaagtag