



Reference: PTXBC004549

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DUS3; tRNA-dihydrouridine(47) synthase [NAD(P)(+)]-like; dihydrouridine synthase 3 like		

For research use only.

GENE DESCRIPTION	dihydrouridine synthase 3-like (<i>S. cerevisiae</i>)
GENE ID	56931
ORIGIN SPECIES	Human
ACCESSIONS	BC004549
VECTOR	pDONR223
NCBI SYMBOL	DUS3L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1