

C14orf129-chromosome 14 open reading frame 129 Gene

Reference: PTXBC004818

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf129; HSPC210; GSK3-beta interaction protein; GSK3B interacting protein		
For research use only.			

GENE DESCRIPTION	chromosome 14 open reading frame 129
GENE ID	51527
ORIGIN SPECIES	Human
ACCESSIONS	BC004818
VECTOR	pDONR223
NCBI SYMBOL	C14orf129
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaacagactgtaatccatggagctaaagcagatgtcaggtatgaagaaggttcagagctgaacggtttgaaggaactgacatgaaagacatgaggctgaagctgaagcagttgtaaatgatgttctctttgctgtaacaacatgtttgtctcgaaagcctgcggtgtcggatga
tggccttatatcaatgtggaacaaaggaaagaaacagatatgcttagaactcactgaagcagggtcaaggtaggtgctatgctttgaccaggtagatgatcatttacagactccctaccatgaacagctactcctgtgttgatactcagcccgctaccgagaagcatttga
aacgcactgctcaagactggaagctttgaaaagagatggacagtcatga