

PAGE1-P antigen family, member 1 (prostate associated) Gene

Reference: PTXBC004861

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AL5; CT16.3; GAGE-9; GAGEB1; PAGE-1; P antigen family member 1; G antigen 9; G antigen family B member 1; G antigen, family B, 1 (prostate associated); P antigen family, member 1 (prostate associated); prostate associated gene 1; prostate-associated gene		

For research use only.

GENE DESCRIPTION	P antigen family, member 1 (prostate associated)
GENE ID	8712
ORIGIN SPECIES	Human
ACCESSIONS	BC004861
VECTOR	pDONR223
NCBI SYMBOL	PAGE1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggttttctaagaagattaatctatcggcgtagaccaatgatctatgtagaatctctgaggagtcagtgatgagcaacctgacgaagtggaaacaccaactcaagtcaggattctacacctgctgaagagagagaggatgaggagcatctgcagctcaaggcaggagcctgaagctgatagccagggaactggttcagccaaagactgggtgtagcttgagatggtctgatacgaagggtgtgctgcgaaatgaagagcagatgaaactgcccgcaagaaggccagagcctgaagcggatagccagggaacagggttcaccgaagactgggtgtgagcgcggagatggtcctgatgtccaggagttggcctgccaaatccagaggaggtgaaaacacctgaggaagatgaagggaatcacagccttaa