



Reference: PTXBC004898

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYP-33; CYP33; peptidyl-prolyl cis-trans isomerase E; PPIase E; cyclophilin-33; peptidylprolyl isomerase E (cyclophilin E); rotamase E; peptidylprolyl isomerase E		

For research use only.

GENE DESCRIPTION	peptidylprolyl isomerase E (cyclophilin E)
GENE ID	10450
ORIGIN SPECIES	Human
ACCESSIONS	BC004898
VECTOR	pDONR223
NCBI SYMBOL	PPIE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]