

RPL36-ribosomal protein L36 Gene

Reference: PTXBC004971

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L36; 60S ribosomal protein L36; ribosomal protein L36		
For research use only.			

GENE DESCRIPTION	ribosomal protein L36
GENE ID	25873
ORIGIN SPECIES	Human
ACCESSIONS	BC004971
VECTOR	pDONR223
NCBI SYMBOL	RPL36
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctacgtaccctatggccgtggcctcaacaagggccacaagtgaacgaagtgagcaagccagccagccgagccggtgtgacaaacacacaaagtcgtgcgggacatgattcgggaggtgtgtggtttgccccgtacgagcgcgccatgga
gttactgaaggtctccaaggacaaacggcctcaaatatcaagaaaagggtgggacgcacatccgcgaaggaagcgggagagctgagcaacgtactggcccatgaggaaagctgctgccaagaaagactga