



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

SLC25A1-solute carrier family 25 (mitochondrial carrier; citrate transporter), member 1
Gene

Reference: PTXBC004980

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CTP; D2L2AD; SEA; SLC20A3; tricarboxylate transport protein, mitochondrial; citrate transport protein; solute carrier family 20 (mitochondrial citrate transporter), member 3; solute carrier family 25 (mitochondrial carrier; citrate transporter), member 1;		

For research use only.

GENE DESCRIPTION	solute carrier family 25 (mitochondrial carrier; citrate transporter), member 1
GENE ID	6576
ORIGIN SPECIES	Human
ACCESSIONS	BC004980
VECTOR	pDONR223
NCBI SYMBOL	SLC25A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1