



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MRLC2-myosin regulatory light chain MRLC2 Gene

Reference: PTXBC004994

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	myosin regulatory light chain MRLC2; MRLC2; MLC-B; myosin regulatory light chain 12B; MLC-2; MLC-2A; MLC20; SHUJUN-1; myosin regulatory light chain 2; myosin regulatory light chain 2-B, smooth muscle isoform; myosin regulatory light chain 20 kDa; myosin,		

For research use only.

GENE DESCRIPTION	myosin regulatory light chain MRLC2
GENE ID	103910
ORIGIN SPECIES	Human
ACCESSIONS	BC004994
VECTOR	pDONR223
NCBI SYMBOL	MRLC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr