



Reference: PTXBC005000

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CFIM; CFIM68; HPBR11-4; HPBR11-7; cleavage and polyadenylation specificity factor subunit 6; CPSF 68 kDa subunit; cleavage and polyadenylation specific factor 6, 68kDa; cleavage and polyadenylation specificity factor 68 kDa subunit; cleavage factor Im com		

GENE DESCRIPTION	cleavage and polyadenylation specific factor 6, 68kDa
GENE ID	11052
ORIGIN SPECIES	Human
ACCESSIONS	BC005000
VECTOR	pDONR223
NCBI SYMBOL	CPSF6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1