



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MRPL11-mitochondrial ribosomal protein L11 Gene

Reference: PTXBC005002

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-113; L11MT; MRP-L11; 39S ribosomal protein L11, mitochondrial; mitochondrial ribosomal protein L11		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L11
GENE ID	65003
ORIGIN SPECIES	Human
ACCESSIONS	BC005002
VECTOR	pDONR223
NCBI SYMBOL	MRPL11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaaaagctcggcggggcccccggggcctcaggaaagcccgaggtctggcgggtgtaactggggcgcagctctggcggcagctccggggccccactagggccagctctgggtcagagaggcgctttccatcaacagattttgcaagggattcaatgagaggacaaagg
acatacaagggaaggcatctcttgcctaccaagattttagtggaagcctgacagacatttgaataaagattggagcagcccatctgtttctactctctgaaggcagcagctggggatgaaagggggccggcgcaaacaggggaaagagctggcaggcctgggtacccttgaagcatgtgtatga
gattgccgcctacaaagctcaggatgagggcatttggccctcagagatgtaccctgtctgttgcctccatcctggctgtctctctctcctcgcgtgtgtaggagcctcagttcagaagaagcttcagcgtttccagaaggaacgagccatctctctgtcgtcctagaagag

More info: Please contact us for any question: contact@proteogenix.fr