



Reference: PTXBC005005

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C17orf80; HLC-8; MIG3; SPEP1; cell migration-inducing gene 3 protein; human lung cancer oncogene 8 protein; lung cancer-related protein 8; sperm-expressed protein 1; chromosome 17 open reading frame 80		

For research use only.

GENE DESCRIPTION	chromosome 17 open reading frame 80
GENE ID	55028
ORIGIN SPECIES	Human
ACCESSIONS	BC005005
VECTOR	pDONR223
NCBI SYMBOL	C17orf80
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr